W NORELCO COOL SKIN



Ultra Close, Ultra Comfortable... Guaranteed!

60-Day Money-Back Guarantee

PHILIPS



7735X

Rechargeable Cordless Razor

FEATURES

- MOISTURIZING SHAVING SYSTEM The Norelco Cool Skin Razor redefines the way you shave. Its unique combination of Improved Lift and Cut® technology with NIVEA for Men shaving lotion guarantees a closer, smoother, more refreshing shave with less irritation
- ULTRA CLOSE NIVEA for Men shaving lotion prepares and softens the beard hair so that the Improved Lift and Cut® dual blade system can cut closer
- ULTRA COMFORTABLE Unique Glide Rings are specially designed to evenly spread and hold lotion in contact with skin longer for smoother gliding and less irritation. NIVEA for Men shaving lotion softens the skin and hair so that razor heads can glide easily over the face for a smoother, more comfortable shave
- **ULTRA REFRESHING** Shave in or out of the shower. NIVEA for Men shaving lotion and gel leave your face feeling smooth and refreshed. Lotion and gel are vitamin enriched to moisturize, and the alcohol-free formulas will not dry or irritate the skin
- Improved Lift and Cut® Technology dual blade system gently lifts hair to cut below skin level
- Patented Reflex Action® adjusts to every curve on your face for flexible, smoother shaving
- Individually Floating Heads align the razor-sharp blades closer to your skin for exceptional closeness
- Stainless Steel Self-Sharpening Blades extend the life of the razor heads
- Quick Clean System Raised head design allows for easy flow through cleaning under
- Full-Width Trimmer defines and shapes moustache and sideburns
- Water Resistant Shave in or out of the shower!
- Shaving Lotion Applicator Use just the amount of lotion you like
- POWER SOURCE
- Rechargeable Cordless for Cordless Operation
- 8-Hour Full Charge for up to 40 minutes of cordless shave time
- LED Battery Charge and Low Indicator
- Automatic Worldwide Voltage

OTHER CONVENIENT FEATURES

- Anti Slip-Grip, Hair Collection Chamber and Electronic On/Off Button
- Charging Stand, Deluxe Travel Pouch, Protective Razor Head Cap, Shower Cartridge Caddy and Hanging Option for convenient in shower use
- 2 Nivea for Men Shaving Lotion and 2 Nivea for Men Gel Cartridges Included
- Norelco 60-Day Money-Back Guarantee
- Full Two Year Warranty

REPLACEMENT ACCESSORIES (sold separately)

HQ177 - Cool Skin 7000 Series Razor Replacement Heads

HQ170 - Replacement Nivea for Men Shaving Lotion Cartridges (5 cartridges per pkg.)

HQ171 - Replacement Nivea for Men Shaving Gel Cartridges (5 cartridges per pkg.)

TECHNICAL DATA

UPC Code Number: 0 75020 19773 3

Dress Carton Dimensions (in.): * 6.69 x 4.73 x 10.55

Cube (ft.): * .19

Weight (lb.): * 1.77 lb. Master Carton Quantity: 3

Weight (lb.): * 5.31 lb.

Master Carton Dimensions (in.): * 7.80 x 13.07 x 11.18

Cube (ft.): * .66

Norelco Consumer Products Company A Division of Philips Electronics North America Corporation Norelco, Reflex Action and Lift and Cut are Registered Trademarks of Philips Electronics North America Corporation. Cool Skin, Quick Clean System and Glide Ring are Trademarks of Philips Electronics North America Corporation. PHILIPS is a Registered Trademark of Koninklijke Philips Electronics N.V.















