

PHILIPS

Saeco

Maintenance
Accessories

Milk Circuit Cleaner

CA6705/99



Perfectly cleans the milk circuit

Keeps milk tubes hygienic and clean

With Saeco milk circuit cleaner you can sanitize the milk circuits of your espresso machine or milk frother. It ensures effective removal of milk residue from all components. We recommend to carry out this cycle once per month.

Effective protection of your espresso appliance

- Protects the milk circuits from milk residue clogging

Prolong lifetime

- Regular cleaning prolongs lifespan of milk frothing parts

Sensational milk foam quality

- Improve the taste of your coffee specialty

Highlights

Protects milk circuits

Protects the milk circuits from milk residue clogging

Prolong lifetime

Regular cleaning prolongs the life of all milk frothing parts (e.g. tubes, Cappuccinatore, Panarello).

Perfect coffee taste



Thanks to the milk circuit cleaner cleaning solution you can effectively remove milk residues that might affect milk taste.

Specifications

Packaging

Quantity: 10 bags of 2 gr.

Weight and dimensions

Weight of product: 0.1 kg

