

PHILIPS

Garment Steamer

ClearTouch

2000 W

Hang&lock, 3 levels

Hanger

Pleat maker, glove



GC536/35



Revive your delicate clothes with steam

Designed for powerful performance and convenience

The new Philips ClearTouch garment steamer is designed for safe and easy crease removal on delicate & difficult to iron clothes. Its unique design, smart features and powerful steam make refreshment of your clothes safe and easy.

Quick crease removal

- 3 steam levels
- Extra powerful steam
- XL Stainless steel steam plate for faster results

Easy to use

- Detachable, transparent water tank with hygienic water inlet
- Easy Rinse descaling function
- Hanger with hang&lock feature for stability during steaming

Maximum safety

- Safe to use on delicate fabrics, like silk
- Glove for extra protection during steaming

Refreshes clothes without washing or dry cleaning

- Removes cigarette, food and body odors
- Steam kills up to 99.9% of bacteria*

Highlights

Extra powerful steam



Powerful steam is blown through the nozzles, enabling you to remove creases with only a few strokes.

3 steam levels



Set your preferred steam level for optimal results on different kinds of garments.

XL Stainless steel steam plate



The steamer head has an XL stainless steel steam plate which delivers faster results.

Safe on all garments



The steamer is safe to use on all fabrics. It is a great solution for delicate fabrics like silk.

Detachable water tank



1.2 L detachable water tank with special water inlet for extra hygiene.

Hanger with hang&lock



Unique hang&lock feature locks the garment hanger and enables stability for more convenient steaming. The removable hanger comes with recesses for dresses and trouser clips.

Easy Rinse descaling function



Prolong the lifetime of your appliance by using the Easy Rinse descaling function regularly.

Steam kills bacteria



Hot steam kills up to 99.9% of bacteria on clothes and helps delaying washing or dry cleaning*.

Removes unwanted odors



Powerful steam removes cigarette, food and body odors.

Glove for extra protection



Included is a glove that protects your hand during steaming.

Specifications

Easy to use

Power cord length: 1.6 m
Water tank capacity: 1200 ml
Detachable water tank
Refill any time during use
Safe on all ironable fabrics: Even delicates like silk
Silicone steam hose
Soleplate name: XL steam plate
Special water inlet: For extra hygiene
Tap water suitable

Accessories included

Adjustable pole
Glove for extra protection
Pleat maker

Fast crease removal

Continuous steam: 40 g/min
Power: 2000 W
Voltage: 220-240 V
Variable steam: 3 levels

Green efficiency

User manual: 100% recycled paper

Scale management

Descaling and cleaning: Easy Rinse

Size and weight

Packaging dimensions (WxHxL):
40x45x33.3 cm
Product dimensions (WxHxL): 33x182x33 cm
Total weight with packaging: 5.7 kg
Weight of iron: 3.7 kg

