



Philips
Steamer

9 L 900 W

Manual timer



HD9120/55

Discover fuller flavours

Add your own taste with the flavour booster

This steamer with unique Flavour Booster improves food with the delicious aroma of herbs and spices. Smart features such as the 60-minute timer, external water inlet and dishwasher-safe parts all make steaming as easy as can be.

Discover fuller flavours

- Flavour Booster improves taste with delicious herbs and spices

Only a little effort

- 60-minute timer
- Dishwasher-safe parts
- External water inlet

Smart features

- Cord storage
- Removable 500 ml drip tray is easy to pour from and clean
- Stackable steaming tiers

PHILIPS

sense and simplicity

Highlights

Flavour Booster



The unique Flavour Booster adds a delicious aroma of herbs and spices, bringing even more taste to steaming. Simply pop your favourite herbs and spices into the booster, and let steam do the rest. The heat from the steam releases delicate aromas from the herbs and spices, which thoroughly infuse the food with their mouth-watering flavours.

60-minute timer



60-minute timer with ready signal and auto shutoff.

Dishwasher-safe parts



Dishwasher-safe parts make cleaning easy.

External water inlet



External water inlet for refilling reservoir during use.

Stackable steaming tiers



Stackable steaming tiers use less space for storing.

Specifications

General specifications

- Number of tiers: 3 pcs
- Water level indicator
- Overheat and dry boil protection
- Cool-touch handgrips
- Non-slip feet
- Power-on light

Technical specifications

- Power: 900 W
- Voltage: 220-240 V
- Frequency: 50/60 Hz
- Cord length: 100 cm
- Capacity tiers: 2.5/2.6/3.5 l
- Capacity water tank: 1.1 l

Design specifications

- Material: Plastic
- Colour(s): White with orange accents

Accessories included

- Rice/soup/food container: 1.2 l
- Egg rack

Weight and dimensions

- Weight of the appliance: 2.3 kg
- Product dimensions (W x H x D): 312 x 450 x 228 mm

