

# PHILIPS

## Digital scale



HF8000



## Give your body an accurate voice

### Fast response

Philips electronic scale gives you an accurate weight measurement. It handles weights up to 150 kg with a precision of 100 g. Stylish, convenient and easy to use.

#### **Professionally accurate measurement**

- Measures up to 150 kg with 100 gram precision
- 4-sensor technology for reliable and accurate measurement
- Adjusts for gravity differences around the world

#### **Fast response**

- Measurement results within seconds
- Instantly ready to use: simply step on

# Specifications

## Features

Easy to clean  
Immediately ready to use  
Anti-slip safety

## Technical specifications

Maximum weight: 150 kg  
Minimum weight: 10 kg  
Resolution: 100 g  
Unit of measure: kg / st / lb  
Gravity compensation: 6 zones  
Technology: 4-sensor  
Display: LED (can also be used in dark places)

Color(s): white

Dimensions: 326 x 326 x 30 mm

Batteries: 4 x LR6 AA 1.5V

## Accessories

Manual

