



Philips  
Digital scale



HF8005

# Give your body an accurate voice

## Weight management

Philips digital scale gives you accurate weight measurement up to 180 kg. Its advanced weight management features make it easy to keep track of changes in your weight and help you on your way to good general health.

### **Professionally accurate measurement**

- Measures up to 180 kg with 100 gram precision
- 4-sensor technology for reliable and accurate measurement
- Adjusts for gravity differences around the world

### **Weight management**

- Shows your weight history over time
- Assists in working towards and maintaining your ideal weight
- Supports multiple users

### **Easy readability**

- Large and high contrast display for easy readability

### **Durable, easy care design**

- Durable and easy to clean stainless steel design

**PHILIPS**  
sense and simplicity

# Specifications

## Features

- Weight management
- Weight history
- Multiple users
- Easy-to-read display
- Durable design
- Easy to clean
- Immediately ready to use

## Technical specifications

- Maximum weight: 180 kg
- Minimum weight: 10 kg
- Resolution: 100 g

- Unit of measure
- Gravity compensation: 6 zones
- Technology: 4-sensor
- Display: LCD
- Color(s): Grey (brushed stainless steel)
- Material: Stainless steel
- Dimensions: 298x298x40 mm
- Batteries: 4 x LR6 AA 1.5V

## Accessories

- Batteries: included
- Carpet feet: 4 pieces
- Manual



Issue date 2010-12-08

Version: 2.1.1

12 NC: 0000 000 00000

© 2010 Koninklijke Philips Electronics N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips  
Electronics N.V. or their respective owners.

[www.philips.com](http://www.philips.com)