

Philips Viva Collection
Juice & Co

300 W

Juicer, blender & filter



HR1842/55



Homemade drinks instantly

Create healthy smoothies and juices

300 W blender and juice extractor with 2 speed settings and pulse, QuickClean function, detachable knife unit, fruit filter and 1.5 l blender jar.

Makes more juice

- Makes more juice thanks to micro mesh filter

Prepare fruit juices and smoothies

- Juicer and blender all in one
- 2 speeds and pulse

Clear juice from your blender

- Filter for clear juice

PHILIPS
sense and simplicity

Specifications

Accessories

- Blender jar
- Cleaning brush
- Filter
- Juice jug

Design specifications

- Color(s): Bright white with orange accents
- Material housing and clamps: PP
- Material lid: SAN
- Material pusher: PP

General specifications

- Integrated cord storage
- Safety clamps
- Speeds: 2
- Suction feet

Technical specifications

- Cord length: 1,2 m
- Frequency: 50/60 Hz
- Power: 300 W
- Voltage: 220-240 V

Highlights

Juicer and blender all in one



You can prepare any juice or smoothie easily.

Multiple speed control



Fast processing and cleaning.

Filter for clear juice



To make soy bean milk and fruit juices in the blender.

