

PHILIPS

WALITA

Steamer

Daily Collection

5 L

900 W

Aroma Infuser



RI9116/01

Pure taste of steamed food

With Aroma infuser

Steam with perfect settings for every meal.

A healthy way to cook tasty meals

- Aroma Infuser adds more taste with delicious herb & spices
- Healthy steaming keeps the nutrition in the food
- Food set timer for fish, vegetables, rice and more

Easily cook a complete meal in one go

- 2 baskets to steam several ingredients at once
- 5L capacity
- Egg rack
- Steaming bowl for soup, stew, rice and more

Smart and convenient handling

- Anti drip lid
- Parts are dishwasher safe
- Easy to fill with water
- Stackable baskets for easy storage

Highlights

2 baskets



2 baskets to steam different ingredients at the same time.

Easy waterfill



Easily refill the reservoir while in use.

Stackable baskets



Stackable baskets for easy storage

5L capacity



5L capacity

Egg rack



Fits 6 eggs per basket.

Steaming bowl



Steaming bowl for soup, stew, rice and more

Anti drip lid



Anti drip lid

Healthy steaming



Healthy steaming keeps the nutrition in the food

Aroma Infuser



The unique Aroma Infuser of the Philips steamer adds delicious aroma of herbs and spices, bringing even more taste to steaming. Simply pop your favourite herbs and spices into the booster, and let steam do the rest. The heat from the steam releases delicate aromas from the herbs and spices, which thoroughly infuse the food with their mouth-watering flavours.

Dishwasher safe parts



Parts are dishwasher safe

Food set timer



Food set timer for fish, vegetables, rice and more

Specifications

Design specifications

Material: Plastic
Color(s): Black

Accessories included

Egg rack
Rice/soup/food bowl

General specifications

Non-slip feet
Power-on light
Cool-touch handgrips
Water level indicator
Number of tiers: 2 pcs
Overheat and dry boil protection
Cord storage

Technical specifications

Power: 900 W
Voltage: 220-240 V
Frequency: 50/60 Hz
Cord length: 80 cm
Capacity (max): 2.6 + 2.9 L
Capacity water tank: 1 L

