

PHILIPS

AVENT

Baby gift set

Infant Starter Set



SCD265//01



Helps settle your baby, especially at night* Less colic**

The Philips Avent Infant Starter Set is the ideal set to get started with bottle feeding. The set includes 4 Classic feeding bottles (2X 40z, 2X 90z), bottle and nipple brush, milk powder dispenser, translucent pacifier and trainer handles.

Vents air away from your baby's tummy

- Clinically proven anti-colic system

More comfortable feed at baby's chosen pace

- Easy latch on due to the unique valve on the nipple

BPA-Free

- Made from Polyether Sulfone (PES)

Benefits and features:

- Bottles clinically proven to reduce colic
- Curved brush head and molded tip
- Holds enough powdered formula for three 260 ml/ 9 oz feeds
- Orthodontic pacifier

Highlights

Bottle and nipple brush



Specially curved brush head and molded handle-tip reach the corners of all types of bottles, nipples and feeding equipment for thorough cleaning.

Bottles clinically proven



Unique anti-colic system consists of the Philips Avent nipple and adaptor ring allowing air to flow into the bottle instead of your baby's tummy for better digestion.*

Easy latch on



Unique valve on the nipple flexes to your baby's feeding rhythm. Milk will only flow at the pace chosen by your baby to minimize overeating and spit-up, burping and gas

Holds enough formula for three 260 ml/ 9 oz feeds



The Philips Avent milk powder dispenser holds 3 pre-measured doses of milk powder - ideal for travel

Orthodontic pacifier



The orthodontic, collapsible nipple of the translucent pacifier respects the natural development of baby's palate, teeth and gums. All our pacifiers are made from silicone and are taste and odour-free.

Unique anti-colic system



As your baby feeds, the unique valve on the nipple flexes to allow air into the bottle instead of your baby's tummy.*

This product is BPA-Free



BPA is Bisphenol-A. This product uses Polyether Sulfone (PES), which is BPA-Free and has a natural honey-colored tint.

Specifications

Material

BPA free*

What is included

Bottle and nipple brush: 1 pcs

Milk powder dispenser: 1 pcs

Baby bottle: 4 pcs

Translucent pacifier: 1 pcs

Trainer cup: 1 pcs

Bottle

Material: BPA free

Design

Bottle design: Wide neck, Ergonomic shape

Ease of use

Bottle use: Easy to clean, Easy to hold

Development stages

Stage: 0-6 months

Functions

Anti-colic valve: Two piece anti-colic system

Latch on: Easy latch on, Promotes nipple acceptance

Nipple: Flexes to feeding rhythm

