



42 codes:	421335400011
Product:	LBL Teat Instructions P&R DB
Languages:	IL (1 lang)
To print colours:	cyan magenta yellow PMS647C PMS3005C
Profile/dimensions:	W30mm x H45mm (keyline)
Material	White
Same profile as:	421335400030
Replaces item no:	421335400010
Point 6 ref:	9586-TeatInstrIL_0011
Date:	7th September 2007
Created by Point 6 Design. +44 (0)20 8962 5880	

P6 PROOF DETAILS

Teat Instruction Label – IL