

# PHILIPS

## AVENT

### Teat

Classic Bottles

Classic

Thick Feed

6m+



SCF636/27



## Good feeding, good sleeping

Helps settle your baby, especially at night\*\*

With the Philips AVENT Classic teat, your baby controls milk flow, which can help reduce overeating and spit-up. The unique skirt flexes to allow air into the bottle and promotes healthy, active feeding and reduces fussing\*\* and colic.\*\*\*

#### More comfortable feed at baby's chosen pace

- Easy latch on due to the unique skirt on the teat

#### Reduce fussing and discomfort

- Clinically proven to significantly reduce fussing\*\*

#### Air goes into the bottle, not baby's tummy

- Unique two-piece anti-colic system

#### BPA Free (0% Bisphenol A)

- BPA free\* teat

#### Other benefits

- Different flow rate teats available
- Compatible with Classic bottle


# Highlights

## Easy latch on



Unique skirt on the teat flexes to your baby's feeding rhythm. Milk will only flow at the pace chosen by your baby for more comfort.

## Helps settle your baby



Sleep and nutrition are vital to your baby's health and happiness. A randomised clinical trial was carried out to see whether infant feeding bottle design affects "infant behaviour". The Philips AVENT Classic baby bottle was shown to significantly reduce fussing by approximately 28 minutes a day as compared to the comparator bottle (46 min vs. 74 min,  $p=0.05$ ). This was especially true at night.

## Unique anti-colic valve



The unique two-piece anti-colic system consists of the Philips AVENT teat and adapter ring. As your baby feeds, the unique skirt on the teat flexes to allow air into the bottle instead of your baby's tummy.

## BPA free\*



This teat is made from silicone – a BPA-free\* material

## Fit your baby's growing needs



The Philips AVENT Classic bottle offers different flow rates to keep up with your baby's growth. Remember that age indications are approximate as babies develop at different rates. All teats are available in twin packs.

## Compatible with Classic bottle

We advise to use Philips AVENT Classic teat only with the Philips AVENT Classic bottle.

# Specifications

## Country of origin

England

## Material

Teat: Silicone, BPA free\*

## What is included

Soft, Thick Feed Teat: 2 pcs

## Design

Teat design: Unique skirt, Easy latch on

## Teat

Flow speed: Fast flow

Holes: Slot hole

Months: 6m+

## Development stages

Stages: 0-6 months

## Functions

Latch on: Easy latch on, Promotes teat acceptance

Teat: Flexes to feeding rhythm

Anti-colic valve: Two piece anti-colic system



\* Following EU directive (2011/8/EU)

\*\* A clinical study showed that at two weeks of age, babies showed less fussing than babies fed with another bottle. (Study conducted by the Institute of Child Health, London. 2008.)

\*\*\* A clinical study showed that at 2 weeks of age, infants fed with an AVENT bottle showed less fussing than infants fed with another leading bottle, especially at night.