

PHILIPS

AVENT

Natural teat

2 pieces

Medium flow

3m+

SCF653/27



The natural way to bottle feed

Avent teat with petal design

Our new teat helps to make bottle feeding more natural for your baby and you. The teat features an innovative petal design for natural latch on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Easy to combine breast and bottle feeding

- Natural latch on due to the wide breast-shaped teat

A more comfortable and contented feed for baby

- Unique Petals for a soft, flexible teat without collapse

Vents air away from your baby's tummy

- Advanced anti-colic system with innovative twin valve

BPA Free (0% Bisphenol A)

- BPA-free teat

Other benefits

- Different flow rates for the most comfortable feed
- Compatible with Philips Avent Natural feeding bottle

Highlights

Natural latch on



The wide, breast-shaped teat promotes natural latch-on similar to the breast, making it easy for your baby to combine breast and bottle feeding.

Unique comfort petals



Petals inside the teat increase softness and flexibility without teat collapse. Your baby will enjoy a more comfortable and contented feed.

Advanced anti-colic valve



Innovative twin valve designed to reduce colic and discomfort by venting air into the bottle and not baby's tummy.

BPA free



This teat is made from silicone - a BPA-free material (Following EU directive 2011/8/EU)

Different flow rates

Philips Avent offers four different flow rates to keep up with your baby's growth. Remember that age indications are approximate as babies develop at different rates. Philips Avent offers the newborn flow teat with 1 hole (0M+), the slow flow teat with 2 holes (1M+), the medium flow teat with 3 holes (3M+) and the fast flow teat with 4 holes (6M+). All teats are available in twin packs.

Compatible with Natural bottle

We advise using the Natural feeding teats with Natural bottles only.

Specifications

Material

Teat: Silicone, BPA-free*

What is included

Soft, Medium Flow teat: 2 pcs

Design

Teat design: Breast-shaped teat, Unique comfort petals, Extra wide

Teat

Flow speed: Medium flow

Holes: 3 holes

Months: 3m+

Development stages

Stage: 0-6 months

Functions

Teat: Unique comfort petals, Extra-soft and flexible teat

Latch on: Easily combine breast and bottle, Natural latch on

Anti-colic valve: Advanced anti-colic system

