

Philips  
Multimedia Speakers 2.1



SPA1315



## Let us entertain you

With 20 W sound output and dynamic bass reproduction, these speakers are perfect for playing music, games, movies and online videos on both MAC and PC. They also come with a handy wired volume control.

### Listen with pleasure

- Dynamic bass
- 20 W sound output power
- Convenient wired volume and power control

### We can play with any media

- We work with MAC and PC
- Enjoy music, games, movies, TV and online videos

**PHILIPS**  
sense and simplicity

# Specifications

## Sound

- Bass enhancement: Bass reflex
- Frequency response: 140–15k Hz Hz
- Music Power: 20 W
- Output power (RMS): Per Satellites: 2 x 2.5 W; Subwoofer 5 W
- Wires: Fixed Wires
- Volume Control: Analogue Volume Control

## Loudspeakers

- Loudspeaker Enhancement: Magnetic Shielded LSB
- Satellite speaker drivers: 2.5"
- Satellite speaker freq range: 200–15k Hz
- Subwoofer driver: 4"
- Subwoofer freq range: 80–240 Hz

## Connectivity

- Cable length: Power cord, audio cable, satellite cable: 1.3 M
- Connector: 3.5 mm stereo
- Audio input for data: Stereo Audio (3.5 mm jack) 1 x
- Power: AC 220 V - 230 V 50 Hz

## Convenience

- Power on indication

## Power

- Power LED indicator: Green

## Packaging content

- Number of satellites: 1
- Subwoofer
- 3.5 mm stereo line cable: Fixed
- Quick installation guide

## Outer Carton


- Gross weight: 11.57 kg
- Outer carton (L x W x H): 48.7 x 33.7 x 45.4 cm
- Net weight: 8.36 kg
- Tare weight: 3.21 kg
- EAN: 87 12581 55535 1
- Number of consumer packages: 4

## Packaging dimensions

- Packaging dimensions (W x H x D): 32.5 x 21.9 x 23.5 cm
- Gross weight: 2.68 kg
- Net weight: 2.09 kg
- Tare weight: 0.59 kg
- EAN: 87 12581 55534 4
- Number of products included: 1
- Packaging type: Carton
- Type of shelf placement: Laying

# Highlights

## Dynamic bass



With 20 W sound output and dynamic bass reproduction.

## Convenient wired control



Easy to control music volume

## Enjoy multimedia experience



The perfect accompaniment for your notebook, so you can play multimedia music and videos with great sound quality.



Issue date 2012-02-27

Version: 1.6.7

12 NC: 8670 000 65882  
EAN: 87 12581 55534 4

© 2012 Koninklijke Philips Electronics N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips  
Electronics N.V. or their respective owners.

[www.philips.com](http://www.philips.com)