

PHILIPS

Straightener

Prestige

Keratin-infused plates

Vibrating plates

Ionic care

230°C professional temperature



HP8371/00



Healthy looking style with extra care

Philips Straightener Prestige provides you with ultimate care while styling with keratin-infused ceramic plates and ionic care. Get optimal results and shiny look with gently vibrating plates.

Beautifully styled hair

- Keratin ceramic plates for smooth gliding and shiny hair
- Ionic care for shiny, frizz-free hair
- Up to 230°C professional high heat for perfect results

Less hair damage

- Gently vibrating plates for optimal styling results
- Gentle temperature setting for touch-ups

Ease of use

- Digital settings with adjustable heat to get smooth results
- Fast heat-up time, ready to use in 30 sec
- Ready to use indicator
- Key-lock function for safe and easy storage
- 2.5m long cord
- Included: heat resistant pouch

Highlights

Keratin ceramic plates



Keratin-infused ceramic plates glide smoothly through your hair, for a shiny look.

Gently vibrating plates



Gently vibrating plates spread hair evenly to enhance optimal styling results for all the hair.

Ionic care



Charged negative ions eliminate static, condition the hair and smooth down the hair cuticles to intensify the hair's shine and glossiness. The result is smooth, frizz-free hair with vibrant shine.

230°C professional temperature



Professional high temperature enables you to change the style of your hair to achieve the look you want.

Fast heat-up time



The straightener has a fast heat-up time, being ready to use in 30 seconds.

Digital temperature settings



Digital display with 12 temperature settings up to 230°C gives you absolute control to adjust temperature to your hair type, for smooth results.

Gentle temperature setting



Refresh your style during the day with gentle temperature setting that prevents overheating and cares for your healthy-looking hair.

Ready to use indicator



This shows when the styler has reached the right temperature, so you can instantly see when you're ready to style.

Key-lock function

The plates can be locked together for safe and easy storage.

2.5m cord



The professional 2.5m long cord reinforces the ease of use anywhere you want.

Included: heat resistant pouch



Included: heat resistant pouch

Specifications

Caring technologies

Ionic Care

Technical specifications

Heating time: 30s

Cord length: 2.5 m

Maximum temperature: 230 °C

Features

Ready for use indicator

Ceramic coating

Material plates: Keratin-infused ceramic plates

Swivel cord

Auto shut-off: after 60 min

Hair type

End result: Straight

Hair length: Long, Medium, Short

Hair thickness: Thick, Thin

Service

2-year worldwide guarantee

