



Philips Viva Collection  
Citrus press

**0.6 L**

with pulp selector



**HR2744/40**

## Enjoy fresh juice in seconds

Citrus press for smooth or pulpy juice

A bright and colorful Philips citrus press HR2744/40 with an adjustable pulp selector on the juicer to adjust the amount of pulp in your juice to your personal taste.

### **You choose smooth or pulpy juice**

- Pulp selector

### **Quick and easy cleaning**

- All parts are dishwasher safe

### **Easy serving**

- 600 ml juice jug

### **Easy storage**

- Compact design
- Cord storage

Citrus press  
0.6 L with pulp selector

# Specifications

## Design specifications

- Color(s): Bright white with blue accents
- Material excl. jug: PP
- Material jug: SAN

## General specifications

- Integrated cord storage
- Non-slip feet
- Pulp selector


## Technical specifications

- Capacity jug: 0,6 l
- Cord length: 1,2 m
- Frequency: 50/60 Hz
- Power: 25 W
- Voltage: 220-240 V

HR2744/40

# Highlights

## Pulp selector



For smooth or pulpy juice.

## Juice jug



600 ml juice jug

## Cord storage



Cord storage

**All parts dishwasher safe**  
All parts are dishwasher safe



Issue date 2013-04-12

Version: 5.0.1

© 2013 Koninklijke Philips Electronics N.V.  
All Rights reserved.

Specifications are subject to change without notice.  
Trademarks are the property of Koninklijke Philips  
Electronics N.V. or their respective owners.

[www.philips.com](http://www.philips.com)